

# BinaxNOW<sup>®</sup>

## Filariasis

The test of choice to help eradicate lymphatic filariasis

120 million people, mostly in developing countries, suffer from this disfiguring disease<sup>1,2</sup>

40 million people are seriously incapacitated by it<sup>1,2</sup>

More than a billion people in over 80 countries are at risk of infection<sup>1,2</sup>

When the result is critical...

...so is the test


[www.binax.com](http://www.binax.com)

# BinaxNOW® Filariasis

Reliably detects circulating antigen to *Wuchereria bancrofti*, the cause of 90% of filariasis infections

## Easy to interpret

- Add sample to the test



- Read at 10 minutes



One line = negative



Two lines = positive



## Ordering Information:

- 620-000 BinaxNOW® Filariasis Test (25 test kit)  
(not available in US or EU)

## Other BinaxNOW Tests:

- BinaxNOW® Influenza A & B
- BinaxNOW® Legionella
- BinaxNOW® RSV
- BinaxNOW® Strep A
- BinaxNOW® S. pneumoniae
- BinaxNOW® Malaria

### References:

1. WHO Fact Sheet No. 102, Revised September 2000.
2. The Global Alliance to Eliminate Lymphatic Filariasis, [www.filariasis.org](http://www.filariasis.org)

## The BinaxNOW® Filariasis test has greatly improved diagnosis and facilitated monitoring of treatment and control programs

- Simple, rapid 'card' test – no need for laboratory facilities
- Uses finger-prick blood taken any time of day, without regard to nocturnicity of parasite
- Detects filarial antigen in whole blood, serum or plasma
- Highly accurate, sensitive and specific
- Easy to use
- Results in 10 minutes

## Distributed by:

Unipath Limited, Bedford MK44 3UP, United Kingdom  
Tel: +44 (0)1234 835000 Fax: +44 (0)1234 835009

